

WWW.IGIA.COM



Sleep Remedy™



SCAN FOR MORE
PRODUCTS



About this Product

Micro-Current Technology helps:

- Stimulate Acupuncture points
- Promote Relaxation
- Release Stress
- Enhance Sleep



The product stimulates the central nervous system of the human body through pulses, so that the body produces melatonin and endorphins, and strengthens the alpha wave of the brain.

Anxiety, depression, insomnia, holding your left hand to divert attention, relieve anxiety Press the metal surface of the product into the palm of your hand and hold it

The hand is the second brain of the human body, the left hand controls the right brain: the right brain is responsible for the human body's emotions, consciousness, and artistry, and the right hand controls the left brain: the left brain is responsible for the body's logic, thinking methods, language skills, reading writing, and computing.

Instructions

Turn on the power switch to choose the frequency that suits your body sensation The low frequency is to soothe and relieve the anxiety frequency. Just place the product in the palm of your left hand The high frequency is the frequency of excitement and encouragement. Just place this product in the palm of your right hand People with anxiety and depression are advised

Instructions for Use

- USB
- Default low frequency mode
- Default high frequency mode



or hold their hands to seep.

to hold each hand for 20 minutes each time, morning, mid-night and night
People with insomnia are advised to hold for 20 minutes before going to bed

Warranty

DISCLAIMER NOTICES

You can register for your free warranty on www.CiA.com. Maintenance cost will be charged for product out of warranty. Damage caused by failure to follow the instructions, human factors, force majeure, etc. is not covered by warranty.

DISCLAIMER NOTICES

DISCLAIMER NOTICES

In the process of the preparation of this manual, we have sought to make the content correct and complete. However, there is no guarantee that there are no errors or omissions in the manual. We reserve the right to change specifications of the products mentioned in this manual at any time without prior declaration. No part of this manual shall be reproduced, transmitted, transferred or stored in a searchable system or translated into any other languages without the written authorization of our company. The company is not responsible for the loss caused by non-standard or wrong application.

Thank you for purchasing this style pie product. Our product's 5 warranty period is 12 months starting from the day of purchase. Please complete the warranty card and use for any future warranty service.

Storing Product

Please try to keep this product in a dry room temperature or in a storage box. Do not overcharge as it may damage battery life. If not used for a long time, please wipe the product and put it into the storage box.

...The red light flashes when charging, the red light is always on when it is full.

the default intensity is 2 when it is turned on.

Buttons "+", "-", and "0": It is the gear to adjust the intensity of EMS, there are 3 intensities, and

frequency mode, and cycle in turn.

Short press the power button to switch the EMS mode: after power on, it will enter the low frequency mode by default, short press the power button to switch to the fast

Long press the power button to turn on, the lights are all on; long press the power button again to turn off, the product will automatically shut down after 20 minutes of use.

Warning

The exterior of the product can be cleaned with alcohol or diluted detergent.

Do not use with implanted medical appliances, such as pacemakers, which may cause equipment failure or trigger a hazard.

may cause equipment failure or trigger a hazard.

Do not use this product if you have skin problems such as wounds, eczema,

Do not use this product in a wet place such as a bathroom.

DO not use this product for anything other than its intended purpose.

Do not disassemble, repair, or modify this product without permission. May

cause equipment failure or injury.

Do not give this product to children under the age of six to avoid personal

·Knfu

Product Name	Remedy Sleep
Material	ABS
Battery Energy	150mah
Input	5V=1A
Working Current	0.1A
Working Wattage	0.4W
Net weight	60G
Size	75x46x25cm
Executive Standard	GB4706.1-2005

Specification

Included:
 1 Remedy Sleep Device
 1 User Manual
 1 Charging Cable