



Items Included




AAA Battery x 2



Travel Pouch

Overview



1. LCD Display Screen
2. Power On & Temperature Reading
3. Temperature Memory Recall
4. Silent Mode
5. Sensor

Instructions

1. Power On

Press the temperature reading button (2) once. You will hear a beep sound indicating the power is on.

2. Temperature Reading

Point the thermometer approximately 1.5" away from the dog and then press the temp read button (2) until you hear a tone sound. For accuracy, aim for the part of the body with the least hair (belly, leg or ears).

3. Silent Mode


Press the silent mode button (4) once to disable sounds. Press the button again to enable sounds.

4. Temperature Memory Recall (Up to 35)


When the thermometer is powered on press the memory button (3) to recall up to 35 previous temperature readings.

Maintenance

Back Cover



AAA Battery x 2



The thermometer's lens is the most sensitive part and needs to be kept clean from dirt or liquids. If the lens becomes dirty, don't touch with your finger or other object, otherwise, it will affect its accuracy. Use an alcohol swab or cotton swab moistened with 95% anhydrous alcohol to clean the thermometer casing and the measuring lens.



FCC Information

This equipment complies with FCC RF radiation exposure limits set forth for an uncontrolled environment.

This device complies with Part 15 of the FCC Rules, operation is subject to the following two conditions:

(1) This device may not cause harmful interference.

(2) this device must accept any interference received, including interference that may cause undesired operation.



Warning: Changes or modifications to this unit not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

• **NOTE:** This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules.

These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications.

However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment on and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help. IC Statement

WARNING/DISPOSAL OF UNIT

Do NOT attempt to change the battery. Doing so will void the warranty and could cause injury and damage to your device. To dispose of the unit, follow the instructions below:

- Never dispose of used batteries/electronics with ordinary solid wastes, since they contain toxic substances.
- Always dispose of used batteries/electronics in accordance with the prevailing community regulations.
- Many electronics retailers now recycle batteries and other electronic components at no charge.

For Product and Warranty Information visit www.MyPetMedics.com or E-Mail Support@mypetmedics.com