15. Do not tamper with this pad in any way. There are no user-serviceable parts. If, for any reason, this pad does not function satisfactorily, return it to authorized service center at address given.
14. Carefully examine inner cover and line cord before each use. Discard the pad if inner covering or line cord shows any sign of deterioration.
13. Do not use the cord as a handle.
12. Never pull this pad by the line cord.
11. Do not sit on or crush pad – avoid sharp folds.
10. Do not use pins or other metallic means to fasten this pad in place.
9. Do not use where oxygen is being administered.
8. Never use pad without removable cover in place.
7. Burns can occur regardless of control setting; check skin under pad frequently.
6. Do not use on areas of insensitive skin.
5. This pad is not to be used on or by any individual with certain disabilities, sleeping or unconscious person, or a person with poor blood circulation unless carefully attended.
4. Do not use on an infant or an animal.
3. Do not use while sleeping. Avoid falling asleep while using the product.

2. ALWAYS UNPLUG THE UNIT FROM THE ELECTRICAL OUTLET IMMEDIATELY AFTER USING.
1. Read all instructions.
- DANGER – To reduce the risk of burns, electric shock and fire, this product must be used in accordance with the following instructions:**

DO NOT DESTROY IMPORTANT SAFETY INSTRUCTIONS

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- therapeutic use of this heating pad.
- **Consult your physician should you have any questions regarding** fracture or persistent pain.
 - undetermined calf pain or joint pain, or if you have poor circulation. Consult a physician before using if you have any kind of fungal infection, swelling, tumors, phlebitis and thrombosis, hemorrhages, inflammations of the skin, by a physician: if you have diabetes, tuberculosis, benign or malignant
 - This heating pad should not be used in the following cases unless directed to OFF when done.
 - **This product does not have an automatic heat-off feature. Slide switch**

SAVE THESE INSTRUCTIONS

24. Do not exceed 30 minutes of treatment at one time.
23. Do not use pad on sensitive skin or when taking medication that can cause sensitive skin.
22. Do not use pad underneath bedcovers and other types of blankets. Do not place covers, such as sheets and blankets, above pad. Only the cloth cover provided with the pad should be used.
21. Do not leave pad unattended, especially if children are present.
20. Do not put this pad in the microwave.
19. Place pad on top and not under the part of the body needing heat.
18. Keep the pad away from heated surfaces.
17. Loop cord loosely when storing; tight wrapping may damage cord and internal parts.
16. Do not use this pad with liniment, salve or ointment preparations that contain heat-producing ingredients, as skin burns could result. Use wherever hot applications are desirable for personal comfort, and wherever recommended by your physician or dentist for the relief of pain.
15. Do not use this pad with liniment, salve or ointment preparations that contain heat-producing ingredients, as skin burns could result. Use wherever hot applications are desirable for personal comfort, and wherever recommended by your physician or dentist for the relief of pain.
14. Do not use this pad with liniment, salve or ointment preparations that contain heat-producing ingredients, as skin burns could result. Use wherever hot applications are desirable for personal comfort, and wherever recommended by your physician or dentist for the relief of pain.
13. Do not use this pad with liniment, salve or ointment preparations that contain heat-producing ingredients, as skin burns could result. Use wherever hot applications are desirable for personal comfort, and wherever recommended by your physician or dentist for the relief of pain.
12. Do not use this pad with liniment, salve or ointment preparations that contain heat-producing ingredients, as skin burns could result. Use wherever hot applications are desirable for personal comfort, and wherever recommended by your physician or dentist for the relief of pain.
11. Do not use this pad with liniment, salve or ointment preparations that contain heat-producing ingredients, as skin burns could result. Use wherever hot applications are desirable for personal comfort, and wherever recommended by your physician or dentist for the relief of pain.
10. Do not use this pad with liniment, salve or ointment preparations that contain heat-producing ingredients, as skin burns could result. Use wherever hot applications are desirable for personal comfort, and wherever recommended by your physician or dentist for the relief of pain.
9. Do not use this pad with liniment, salve or ointment preparations that contain heat-producing ingredients, as skin burns could result. Use wherever hot applications are desirable for personal comfort, and wherever recommended by your physician or dentist for the relief of pain.
8. Do not use this pad with liniment, salve or ointment preparations that contain heat-producing ingredients, as skin burns could result. Use wherever hot applications are desirable for personal comfort, and wherever recommended by your physician or dentist for the relief of pain.
7. Do not use this pad with liniment, salve or ointment preparations that contain heat-producing ingredients, as skin burns could result. Use wherever hot applications are desirable for personal comfort, and wherever recommended by your physician or dentist for the relief of pain.
6. Do not use this pad with liniment, salve or ointment preparations that contain heat-producing ingredients, as skin burns could result. Use wherever hot applications are desirable for personal comfort, and wherever recommended by your physician or dentist for the relief of pain.
5. Do not use this pad with liniment, salve or ointment preparations that contain heat-producing ingredients, as skin burns could result. Use wherever hot applications are desirable for personal comfort, and wherever recommended by your physician or dentist for the relief of pain.
4. Do not use this pad with liniment, salve or ointment preparations that contain heat-producing ingredients, as skin burns could result. Use wherever hot applications are desirable for personal comfort, and wherever recommended by your physician or dentist for the relief of pain.
3. Do not use this pad with liniment, salve or ointment preparations that contain heat-producing ingredients, as skin burns could result. Use wherever hot applications are desirable for personal comfort, and wherever recommended by your physician or dentist for the relief of pain.
2. Do not use this pad with liniment, salve or ointment preparations that contain heat-producing ingredients, as skin burns could result. Use wherever hot applications are desirable for personal comfort, and wherever recommended by your physician or dentist for the relief of pain.
1. Do not use this pad with liniment, salve or ointment preparations that contain heat-producing ingredients, as skin burns could result. Use wherever hot applications are desirable for personal comfort, and wherever recommended by your physician or dentist for the relief of pain.

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IB-14733A

LIMITED FIVE-YEAR WARRANTY (U.S AND CANADA ONLY)

Conair will repair or replace (at our option) your unit free of charge for 60 months from the date of purchase if the appliance is defective in workmanship or materials.

To obtain service under this warranty, return the defective product to the service center listed on back, together with your purchase receipt and \$3.00 for postage and handling. California residents need only provide proof of purchase and should call 1-800-3-CONAIR for shipping instructions. In the absence of a purchase receipt, the warranty period shall be 60 months from the date of manufacture.

ANY IMPLIED WARRANTIES, OBLIGATIONS, OR LIABILITIES, INCLUDING BUT NOT LIMITED TO THE IMPLIED WARRANTY OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE, SHALL BE LIMITED IN DURATION TO THE 60-MONTH DURATION OF THIS WRITTEN, LIMITED WARRANTY. Some states do not allow limitations on how long an implied warranty lasts, so the above limitations may not apply to you.

IN NO EVENT SHALL CONAIR BE LIABLE FOR ANY SPECIAL, INCIDENTAL, OR CONSEQUENTIAL DAMAGES FOR BREACH OF THIS OR ANY OTHER WARRANTY, EXPRESS OR IMPLIED, WHATSOEVER. Some states do not allow the exclusion or limitation of special, incidental, or consequential damages, so the above limitation may not apply to you.

This warranty gives you specific legal rights, and you may also have other rights, which vary from state to state.

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Visit us on the web at www.conair.com
©2018 Conair Corporation
East Windsor, NJ 08520 / Glendale, AZ 85307

HP08F
18PA058165

IB-14733A

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©2019 Conair Corporation