



dpl[®]

Oral Care

Light Therapy System

INCLUDED

- Oral Care Control Device
- Interchangeable Mouth Tray
- Power Supply
- Power Cable
- Storage Case
- Silicone Cleaning Cover



WELCOME

Welcome to LED Technologies, Inc and our dpl® Oral Care line. We are thrilled to announce this product as an in-home device to be used at your convenience.

The dpl® Oral Care Light Therapy System uses light emitting diodes (LEDs) in the infrared and red spectrums, which are medically proven to **increase circulation, enhance gum health, and relieve pain**. Non-UV blue lights **increase the effectiveness of whitening gels by removing tannic stains from coffee, tea and wine, and leaving whiter teeth**.

Please read this manual before operating the Oral Care System.

CHARGING INSTRUCTIONS

- **Only use provided Power Supply. Attempts to use a third party power supply could damage device and will void warranty.**
- Charging:
 - Charging typically takes 4-6 hours but may take longer.
 - Plug Power Cable into Power Supply and Oral Care Control Device. Plug Power Supply into power source.
 - Separate the Oral Care Control Device from Inter-changeable Mouth Tray.
 - While charging, red light flashes until fully charged.
 - When fully charged, red light is constant.
- Battery Indicator:
 - When battery is fully charged, green light is on and red light is off.
 - When battery level reaches less than 25% charge, green light begins to flash.
 - When battery level reaches less than 5%, red light begins to flash.
 - When battery level is close to 1%, the LEDs are powered off and the device is powered off.

OPERATING YOUR ORAL CARE SYSTEM

- Push top button to see lights appear.
- Place device in mouth.

TREATMENT

TEETH WHITENING:

- For first week, use once daily once a day for 20 minutes, with teeth whitening gel of choice.
- After first week, use 2-3 times per week or as needed, with teeth whitening gel of choice.

GUM HEALTH:

- Use once daily for a minimum of 10 minutes.

HOW TO CLEAN DEVICE

- Pull apart device and install waterproof Silicone Cleaning Cover over Interchangeable Tray USB connection.
- Wash mouthpiece in lukewarm water then let device air dry.
- Remove waterproof silicon cleaning cover.

CUSTOMER SERVICE

Need further assistance? You can reach LED Technologies, Inc's Customer Service by calling (800) 337-9565, Monday-Friday, 7:30am-5:00pm, Mountain Standard Time, or by emailing customerservice@ledtechnologies.com.

CONNECT WITH US

When you register your **dpl**® Oral Care Light Therapy System device at revivelighttherapy.com, you will be entered into our system to be eligible for exclusive offers, VIP discounts, new product previews, and the opportunity to earn free products through new product trials.

REPLACEMENT PARTS

The following lost or damaged parts can be purchased as replacement items directly from LED Technologies, Inc by calling Customer Service at (800) 337-9565.

- Power supply
- Power cord
- Interchangeable Mouth Tray

DISCLAIMER

This product is not intended to diagnose, cure or prevent disease. LED Technologies, Inc makes no claims, representations or warranties regarding the ability of this product to cure any physical, oral or mental condition. A qualified health professional should be consulted with regard to any condition requiring medical attention.

WARRANTY

LED Technologies, Inc warrants that this product will be free from defects in materials and workmanship for 6 months from the date of purchase by the original owner. This warranty covers the original purchaser only and is not assignable or transferable to another party. Register your warranty online at revivelighttherapy.com or by calling Customer Service at (800) 337-9565.



6000 Greenwood Plaza Blvd Ste 110
Greenwood Village, CO 80111
800.337.9565

Our products are continually
improving. The actual product may
differ from the product shown.