

PROFESSIONAL HAIR REMOVAL DEVICE

Item No. 207824

User Guide



Thank you for purchasing the Professional Hair Removal Device. Please take a few minutes to read this guide and keep it for future reference.

INTRODUCTION

The Professional Hair Removal Device uses intense pulsed light to neutralize hair follicles and break the cycle of hair growth. Each pulse of light targets the melanin deep inside the hair follicle, causing the hair to fall out naturally and prevents new hair from growing back.

FEATURES

- Clinical-grade technology for home use
- Suitable for hairline, cheek, lip, back, armpit, chest, arm, abdomen, bikini line, thigh and calves
- Compact and portable
- Includes razor and safety glasses
- Plugs into an AC outlet

Please note: This device is not suitable for those with dark skin [dark brown to black] and is not effective on white, gray, red or light blonde hair.

DESCRIPTION OF PARTS



HAIR COLOR AND SKIN TONE GUIDE

Please consult the chart below to see if your hair color and skin tone are suitable for treatment with the Professional Hair Removal Device.

Hair Color	Skin Tone					
	White (Caucasian)	Fair (East Asian)	Medium Brown	Light Brown	Brown	Dark Brown/Black
White/Gray	NO	NO	NO	NO	NO	NO
Red	NO	NO	NO	NO	NO	NO
Light Gold	NO	NO	NO	NO	NO	NO
Light Brown	YES	YES	YES	YES	YES	NO
Dark Brown	YES	YES	YES	YES	YES	NO
Black	YES	YES	YES	YES	YES	NO

BASIC OPERATION

IMPORTANT: *Always wear the included safety glasses when operating the Professional Hair Removal Device.*

1. Use an alcohol wipe to clean the flash window. Shave the target area prior to use. Remove any residual hair.
2. Plug the device into an AC outlet. The white indicator light will flash.
3. Press and hold the Power button for 3 seconds. The white indicator light will turn solid and the cooling fan will turn on.
4. Click the Power button to select the desired intensity level [1-5]. The green LEDs will indicate the intensity level.
5. Move the flash window to the target area until the white indicator light starts to flash. Now you are ready to activate the pulse light. Press the Flash button to flash one time. **NOTE: *The white indicator light will flash when in proximity to the skin, due to a built-in sensor in the flash window. If the white indicator light is not flashing, move the device closer to the skin.***
6. Move the device to an adjacent section of the target area and flash again. Continue moving the device and flashing one time until you have covered the entire target area. **NOTE: *After every flash, the device needs to recharge. When the device is ready to flash again, the white indicator light will illuminate [in 1-3 seconds].***
7. Press and hold the Power button to turn the device OFF. Clean the flash window with an alcohol wipe. **NOTE: *When you are finished, cool the target area with ice water, a damp towel or a cooling balm. A toner may also be used to retain moisture.***

TESTING THE DEVICE

Before beginning the hair removal process, test the device by flashing it on the inside of your upper arm 1-2 times, from low to high intensity. Wait 2 hours to see if you have any adverse reactions. If you experience severe reactions [redness, burns, stabbing pain], discontinue use and contact a physician. If slight discomfort occurs, lower the intensity before using again. If the discomfort still occurs, consult a physician.

OPERATING TIPS

- Patience is key to using this device.
- The hair removal cycle will be about 2-3 months.
- Use the device for 1 session every 2 weeks during hair removal cycle.
- Do not use the device multiple times on the target area. Consult the chart below for recommended uses.
- Trim or shave thick hair prior to using the device.
- As the hair growth diminishes, use the device less frequently.

NOTE: *Results may vary.*

RECOMMENDED USES FOR YOUR SKIN TONE

Body Part	Skin Tone					
	White (Caucasian)	Fair (East Asian)	Medium Brown	Light Brown	Brown	Dark Brown/Black
Armpit	Level 3-5 3-4 sessions	Level 3-5 3-4 sessions	Level 3-4 2-3 sessions	Level 2-3 2-3 sessions	Level 2-3 2-3 sessions	Not applicable
Arm	Level 3-5 3-4 sessions	Level 2-4 3-4 sessions	Level 2-4 2-3 sessions	Level 1-3 2-3 sessions	Level 1-3 2-3 sessions	Not applicable
Leg	Level 3-5 3-4 sessions	Level 2-4 3-4 sessions	Level 2-4 2-3 sessions	Level 1-3 2-3 sessions	Level 1-3 2-3 sessions	Not applicable
Around mouth	Level 3-4 3-4 sessions	Level 2-3 3-4 sessions	Level 2-3 2-3 sessions	Level 1-2 2-3 sessions	Level 1-2 2-3 sessions	Not applicable
Face	Level 3-4 3-4 sessions	Level 2-3 3-4 sessions	Level 2-3 2-3 sessions	Level 1-2 2-3 sessions	Level 1-2 2-3 sessions	Not applicable

DEVICE OPERATION NOTES

IMPORTANT: *Always wear the included safety glasses when operating the Professional Hair Removal Device.*

1. Do not flash on or directly near the eyes.
2. Do not clean the device with spray alcohol.
3. Do not wear cosmetics prior to use.
4. Avoid light-sensitive food [i.e., spinach] during operation.
5. For best results, use the device at higher intensity levels. **(Before increasing intensity levels, be sure to perform the skin sensitivity test as described earlier.)**
6. Darker skin tones may not provide immediate skin reactions. Follow the 2-hour waiting period as described above before increasing the intensity.
7. Skin moisturizers and sun block can be used as directed after using the Professional Hair Removal Device.
8. Avoid UV rays [sunlight] during use. Protect skin while in the sun.
9. Do not block the air inlet and the cooling outlet, as this can affect the device's cooling process.
10. Hair stuck to the device may result in blank spots. However, this does not affect the performance of the Professional Hair Removal Device.
11. Not for use by children. Keep out of reach of children at all times.
12. Do not open or attempt to repair this device. There are no user-serviceable parts.

HEALTH WARNINGS

1. **WARNING!** Do not use this product if you are pregnant, under the age of 18, incapacitated or immobile or if you cannot clearly understand the instructions and/or operate the controls.
2. **WARNING!** Do not use this product if you have any heart-related illnesses, suffer from epilepsy, fever, inflammation, dermatitis [acne, eczema, herpes], warts open or recently stitched wounds, serious or progressive illnesses [cancer, AIDS, hepatitis, blood diseases, etc.], blood coagulation problems or mental illness.
3. **WARNING!** Do not use the device on your head, eyes, neck [throat], breasts, genital area or inside your mouth.
4. **WARNING!** Do not store the device where it could encounter vibrations or shock.
5. **WARNING!** Do not store the device in a dusty, salty or sulfuric environment. These elements can have a negative effect on the performance of the Professional Hair Removal Device and decrease life expectancy.

SPECIFICATIONS

Lamp Tube	Xenon quartz tube
Light Source	Intense Pulsed Light
Energy Frequency	1.5-3.9J/cm ²
Focus Size	4.2 cm ²
Wavelength	470nm-1100nm
Lamp Tube Life	600,000 pulses
Adapter	AC100~240V @ 12V DC 3A
Weight	0.5 lbs.
Size	5.5" x 3.3" x 1.9"
Operating Temperature Range	41°F-86°F
Operating Humidity Range	Relative Humidity 30%-60%