

MAGIC TENS OPERATION STEPS

① Install



Connect two buttons on the the electrode pad to the main unit;

② Get ready



Remove the TENS unit from its attached paper, and place it on your neck, shoulder, arm or leg area;

③ Start TENS



Press Power button (indicate the shape or log of the button) for 3 seconds, once beep sound heard, green dot light indicates the unit is ready to work;

④ Intensity Adjustment



Press **+** / **-** Power button to adjust **intensity**, two "beeps " indicate at the lowest / **highest** intensity;

⑤ Mode Adjustment



Single-click on/off button, it will beep once each time to change to a different mode; total 5 modes;

⑥ Finish up



Auto Shut-off: the Magic TENS will shut off automatically in 15mins;
Manually Shut-off: press Power button (show the image) for 3 seconds to shut off.