

AIROSAGE BY FOOT DR.™

Cordless, Portable Leg Massager with Heat Therapy



USER MANUAL

Please read these instructions carefully before using the product

AIROSAGE: Cordless, Portable Leg Massager with Heat Therapy

Product Introduction

AIROSAGE utilizes a specially developed method of massage that uses compression air therapy to revitalize fatigued arms and legs. Upon your first use, you will instantly understand the many benefits of using the AIROSAGE massager that has been designed to help improve circulation. The AIROSAGE relieves minor muscle aches and pains and temporarily helps increase circulation to treated areas.

AIROSAGE helps increase your circulation after a long day on your feet or during a long flight. You can operate it anywhere as it receives its charge from a battery pack. Just secure it around your arm or leg, turn it on and let the compression air massage and soothing heat therapy release tension, stress and pain.



- May help improve circulation.
- May relieve pain in legs or arms due to too much or too little movement
- Sensor regulates air pressure.
- Offers a choice of preprogrammed or manual modes. Both are 15 minutes
- Portable for use anywhere.
- Perfect for frequent flyers, &/or anyone suffering from muscle stiffness.

Instructions

- Start by connecting the controller to the power bank.
- Wrap the device around the area you would like to treat (I.e. Arms or Legs)
- Turn the device on by pressing and holding the “Power” button for 3-5 seconds. A red light will turn on the controller to indicate the device is powered.
- Select the massage option by pressing the “Massage” button. A red light will turn on under “Massage” to indicate this option is powered on
- Select the heat option by pressing the “Heat” button. A red light will turn on under “Heat” to indicate this option is powered on
- To switch off the device, press the “Power” button repeatedly until the red indicator light turns off.



***Please Note:** The device is equipped with a built-in timer and will turn off automatically after 30 minutes.

Product Specification

Main Device

- To turn the power off, remove the batteries from the unit and properly store.
- Always keep clean the unit by using a soft cloth to clean the surface of the unit.
- Shell is non-waterproof. Do not let water enter the product, do not touch with wet hands and don't place or use in moist environment.
- Do not spill liquid on the device
- Do not clean with chemicals

Disclaimer and Warranty

Disclaimer

In the process of the preparation of this manual, we have sought to make the content correct and complete. However, there is no guarantee that there are no errors or omissions in the manual. We reserve the right to change specifications of the products mentioned in this manual at any time without prior declaration. No part of this manual shall be reproduced, transmitted, transcribed or stored in a searchable system or translated into any other languages without the written authorization of our



company. The company is not responsible for the loss caused by non-standard or wrong operation. In the process of the preparation of this manual, we have sought to make the content correct and complete.

Warranty

1. Damage caused by failure to follow the instructions, human factors, force majeure, etc. is not covered by warranty.
2. Maintenance cost will be charged for product out of warranty.
3. You can register for your free warranty on www.IGIA.com.

Warnings

If you are in the care of a physician, consult your physician before using this device.

Safety Precautions:

- Do not use during pregnancy.
- Adult supervision is required when the device is to be used by a handicapped individual.
- Do not use on children, keep out of reach from children.
- Do not use if you have a skin disorder, broken skin, injury, post-surgery or high blood pressure.
- It is always best that you consult with your doctor prior to use of the product if you have medical condition.
- Do not in any way try to alter or service the unit, it will void your warranty. We are not responsible for any misuse other the intended use of the product.
- Do not use in wet conditions or immerse in water.
- Wipe clean with a clean damp cloth.
- Do not use this device over swollen or inflamed areas or skin eruptions or in the presence of unexplained calf pain.
- Do not use this device over the site of or in the vicinity of a broken bone or healing fracture.
- Do not use this device during pregnancy or right after giving birth.
- Do not use this device if you are currently suffering from Deep Vein Thrombosis, cerebral or coronary thrombosis, severe peripheral circulation problems, malignant tumors, heart disorders, severe aneurysms, varices, sensory disorders due to diabetes, lack of sensitivity to temperature, or high fever.
- Consult with your physician before using this device:
 - If you have a pace maker or another implanted medical device.
 - If you are experiencing a menstrual period.



Customer Service:

Phone: 516-653-0666

Email: AIROSAGE1@gmail.com

www.IGIA.com



Join the Club & Save 20% Off Every Order

- As an IGIA™ Savings Club Member, You will receive 20% OFF any item advertised on our website.
- "No Worry" Membership Renewal: A one time fee of \$29.99 for a years worth of savings
- Save BIG on every purchase you make, potentially saving you hundreds of dollars
- We Guarantee our Club will pay for itself, time after time.

Join now & start saving today!

Visit us at **www.IGIA.com**

CONTACT US

Call us 1(516) 653-0666
Email: Sales.Igiany@gmail.com

Register for your FREE warranty

Thank you for shopping with us! We are confident that you will be very satisfied with your purchase. To protect your purchase and register for your warranty, please visit us at www.IGIA.com. Located at the top right hand of the website you will find an icon listed "Warranty Registration"; This is where you will fill in the necessary fields to register and receive your warranty.

We appreciate your support and look forward to your future savings! If you have any questions, please don't hesitate to contact us at: Sales.IgiaNY@Gmail.com

Visit us at www.IGIA.com

0716

LOT#:020720