



SHARPER IMAGE®

FLEX ROLLER
MODULAR MASSAGER

CONGRATULATIONS ON YOUR PURCHASE OF THE FLEX ROLLER MODULAR MASSAGER.

Please read all instructions and warnings prior to use. Keep this manual as it contains important information for future reference.

SHARPER IMAGE® NAME AND LOGOS (INCLUDING THE 77 LOGO AND TOMORROW'S TOMORROW) ARE TRADEMARKS OR REGISTERED TRADEMARKS. DISTRIBUTED BY AND © 2022 MERCHSOURCE, LLC. IRVINE, CA 92618. ALL RIGHTS RESERVED.

FOR CUSTOMER SUPPORT, PLEASE CALL
1-800-374-2744.

PRINTED IN CHINA.


SAFETY WARNINGS

PLEASE READ BEFORE USE.

- Not a children's item. Keep out of the reach of children.
- Do not use this item if you are pregnant or have an injury or illness. Consult your physician if you have a medical condition that may prevent safe use of this item.
- Do not use attachments of any kind other than those that may be included.
- For indoor use only.
- Thoroughly inspect products before every use. If the product appears worn, frayed/splintered, cracked, or broken in any manner, discontinue use and discard immediately.
- If you experience pain while using the massager, stop immediately. Consult a physician before resuming use of the massager.
- Be careful to avoid pinching skin or muscles when assembling and using this massager.
- Keep loose clothing and hair away from the rollers.

WHAT YOU GET



- a. Flex Roller U-shaped handle
- b. T-shaped grips (x2)
- c. Smooth rollers (x4)
- d. Textured rollers (x4)



HOW TO USE ATTACHMENTS


- The Flex Roller is designed to be used in 3 key ways and includes 2 sets of massage rollers.

1. Select between the 2 sets of rollers.
 - Smooth - Soft rollers
 - Textured - Semi-firm rollers
2. Remove T-shaped grip from main U-shaped handle by pressing the orange button at the bottom of the T-shaped grip and sliding it off.
3. Slide roller onto the grey spindle until it is fully seated.
4. To remove roller, slide it off.
5. Reattach T-shaped grip to U-shaped handle by sliding it back until it clicks securely.




MASSAGING ARMS + LEGS

1. Attach both T-shaped grips to the U-shaped handle with the rollers facing inward.
2. Hold each side of the U-shaped handle and flex it to spread the rollers until they fit over your arm or leg.
3. Massage the desired area.



MASSAGING LARGER BODY AREAS

1. Attach both T-shaped grips to the U-shaped handle with the rollers facing inward.
2. Position the rollers flat against the desired area.
3. Move all 4 rollers over the area at the same time.



MASSAGING TARGETED AREAS

1. Remove the T-shaped grips from the handle.
2. Hold the top of the T-shaped grip and use the 2 rollers for a targeted massage.

