





黑白印刷，128G哑粉纸，2折3页，58\*94mm

<div><p>DL004 Solar VBar for Face Operating manual</p></div> <div><p><b>Introduction</b></p><p>Prospera's Solar V bar for is fashionably designed to smooth wrinkles, lift and tighten the skin, shape the face, increase skin metabolism by promoting blood circulation, and reduce the appearance of fine lines. Daily use of the solar V Bar allows your face and skin to look healthier and younger.</p></div>	<div><p><b>Content</b></p><p><b>CONTENTS</b></p><p>Safety Warnings ..... 2</p><p>Product Illustration ..... 2</p><p>Operating the Product ..... 3</p><p>Consumer Warranty ..... 4</p><p>Product Specifications ..... 4</p><p>Warranty Card ..... 5</p></div> <div><p>1</p></div>	<div><p><b>Safety Warnings</b></p><p>Consult your doctor prior to using the product if you have heart, skin and/or other health-related issues.</p><p>1. This product is intended for adult use only.</p><p>2. Do not let babies and children play with the product.</p><p>3. Do not use the product on an infant.</p><p>4. If you experience discomfort, stop using the product immediately.</p><p>5. Do not wash the unit using organic solvent such as thinner to prevent damage.</p></div> <div><p><b>Product Illustration</b></p><p>V-shaped Bar for Close Contact</p><p>WEDM Itched Grooves for Effective Face Lift</p><p>Solar Energy Inducer</p><p>100% Stainless Aluminum Bar</p></div> <div><p>2</p></div>
<div><p><b>Operating the Product</b></p><p><b>Step 1</b> Wash your face, and follow with a small amount of facial serum or moisturizer.</p><p><b>Step 2</b> Roll the VBar in the direction of the arrows illustrated in the above diagram. Repeat each motion 10 times per area.</p><p><b>Step 3</b> Roll the Vbar on both the front and back of your neck in an upward motion. Repeat 10 times for each area. You can use this same technique on other parts of the body, such as the chest, hands and arms, to improve the complexion of the skin in these areas.</p><p>You may feel a minor stinging sensation (activation) as small amounts of micro current generated by the Solar inducer tighten up skin and improve blood circulation.</p><p>Recommend usage: 5-7 times per week</p><p>To order a new health and beauty product, please visit <a href="http://www.prosperacorp.com">www.prosperacorp.com</a>.</p><p>Use coupon code: PROS\$5 for \$5 off any product above \$30.</p></div> <div><p>3</p></div>	<div><p><b>Consumer Warranty</b></p><p>Prospera offers a 90-day standard warranty on massage products. Labor and parts are included during the warranty period. Please provide the model number and dated proof of purchase (sales receipt) when contacting Prospera for warranty service.</p><p>Non-warranty service is provided on a "per incident" basis. The consumer will verify that the product has failed and provide information for servicing the unit. All applicable repairs, parts, shipping, handling, local tax and "per incident" fees will be charged for non-warranty repairs and support calls.</p></div> <div><p><b>Product Specifications</b></p><ul style="list-style-type: none"><li>• 120°V shape roller, 3D contact,</li><li>• Rollers, WEDM itched grooves for effective face lift</li><li>• Direct contact with skin during usage</li><li>• Improve loose skin, eye bags &amp; facial complexion</li><li>• Portable and convenient to use anywhere, anytime</li><li>• Solar energy for easy use</li><li>• Highest grade metal material</li><li>• Current 0.5mA, voltage, 1.5V</li><li>• Power 1W</li><li>• Measurement: 3.5" x1.4" x10.7"</li><li>• Net product weight 2.7 oz</li><li>• Gross Weight: 5.3 oz</li><li>• UPC 850011335024</li></ul></div> <div><p>4</p></div>	<div><p><b>Warranty Card</b></p><div><p>Warranty Registration Form</p><p>Product Name: _____</p><p>Product Model: _____</p><p>Purchase Price: _____</p><p>Date of Purchase: _____</p><p>Retailer: _____</p><p>First Name: _____</p><p>Last Name: _____</p><p>Address: _____</p><p>City: _____ State: _____ Zip: _____</p><p>Phone: _____ Email: _____</p><p>Signature and date: _____</p></div><p>Please mail to:</p><p>Prospera Corporation 392 West Larch Rd Suite 31 Tracy, CA 95304-1644 Tel: (925) 292-5978 Email: <a href="mailto:sales@prosperacorp.com">sales@prosperacorp.com</a></p></div> <div><p>5</p></div>