



# 4-IN-1 BODY TONING FIRMING AND SCULPTOR

Body Sliming

Face Lifting

The combination of 3 technologies makes your body more perfect



**SUCTION + EMS + LED RED LIGHT**

# EMS Function

**Building muscle**

**Promote blood circulation**

**Firming and wrinkle removal**

How it works:

EMS is called electrical muscle stimulation. It stimulates motor neurons for muscle contraction. It is a passive exercise that makes muscles more compact and elastic. At the same time, it can promote blood circulation to a certain extent, so that the production of collagen and elastic fibers can reach a wrinkle firming effect.



Frequency parameter: 10-45Hz

# Our Unique EMS Technology

## Our EMS can stimulate both fast and slow muscles


- Muscles are composed of muscle fibers. Fast and slow muscles need to be trained through EMS to achieve the effect of firming the skin and increasing elasticity.

## Our EMS has both low frequency and medium frequency, and can train two kinds of fibers

- Pulse frequency - different frequencies for different muscle fibers
- Low frequency ( $\approx 20\text{Hz}$ ) stimulates slow-twitch muscle fibers, increasing muscle tolerance and elasticity
- Intermediate frequency ( $\approx 40\text{Hz}$ ) stimulates fast-twitch muscle fibers and increases muscle strength

## We also set the EMS frequency as a step-by-step change, which can enhance the muscle training effect

- Variations in pulse frequency can reduce muscle group fatigue, improve body sensation and training effectiveness




slow-twitch  
muscle fibers

fast-twitch  
muscle fibers

# Our Unique EMS Technology


BCE1



- ◆ **4-fold frequency conversion - combined with mid-low frequency**
  - Effectively trains two types of muscle fibers (red muscle fibers and white muscle fibers)
  - Achieve better muscle training results
  - Comfortable skin nerve adaptation
  - Rich and delicate body feeling

Including low frequency  
and mid frequency  
Incremental frequency  
change

Other Products



- ◆ **General EMS frequency conversion**
  - Trains one type of muscle fiber (red muscle fiber)
  - Easy to cause skin nerve stimulation

Only low frequency  
interleaving frequency  
changes

# Negative Pressure Liposuction Function



Pull up subcutaneous adipose tissue to promote blood circulation, accelerate new skin metabolism, remove cellulite

Negative pressure suction force: 20-40kPa

Promote blood circulation

Relieve body edema

Eliminate cellulite

How it works:

By evacuating the air, the pressure in the suction cup is reduced to absorb the skin tissue, and the microvessels are dilated and congested by physical force, which increases blood circulation, improves metabolism, and finally achieves the effect of eliminating cellulite.



# LED RED LIGHT Function


Promote collagen regeneration

Improve dull skin

Brighten skin colour

How it works:

LED red light is a therapy that achieves cosmetic effects through non-thermal, non-invasive light. It stimulates the mitochondria of cells through the photobioregulation mechanism, increases the proliferation of fibroblasts and the synthesis of collagen, and achieves the cosmetic effect of skin rejuvenation and brightening.



Red wavelength:  $630\pm 5\text{nm}$

# BODY MODE

SUCTION + EMS + RED LIGHT



**Big Suction Cup**  
For legs, stomach, buttocks



**Small Suction Cup**  
for arm

**LIPO function:** 3 levels  
**EMS function:** 3 levels  
**Red light:** Start at power on

According to the acceptable level of the body,  
the level can be adjusted at will, and the suction force and the  
intensity of the EMS can be controlled.  
Comfortably creates a firmer body line.

Mode Key first





# FACE MODE

EMS + RED LIGHT



Round massage brush head design

6 layers of electrodes, comfortable to create a V face

**EMS function:** 3 levels

**Red light:** Start at power on

Adjust the current in multiple gears according to the tolerance of the face

EMS lifts the face while the red light brightens the complexion

Mode Key first



# Excellent minimalist design



Streamlined ergonomic design  
Comfortable to hold


Unique non-slip plaid  
Prevent slipping

Stainless steel massage head  
Safe and durable

Round massage head  
Fits the skin comfortably

Unique non-slip plaid

# Specifications



Product	BCE1 Anti Cellulite EMS & Vacuum Massager
LIPO	Suction power: 50-51 kPa
EMS	Body Mode: 39-115V (3 Levels) Face Mode: 23-45V (3 Levels)
Red Light	Wavelength: 630±5nm
Battery Time	Charging time: 3-4 hours Fully charged usable time: ≥60 minutes
Automatic Shutdown	10 minutes
Charging Method	Type-c
Waterproof	No
Size	Device: 145 * 60 * 69mm Big Suction Cup: 96*63* 33mm Small Suction Cup: 24*24*29mm
Weight	Device: 234g Big Suction Cup: 21g Small Suction Cup: 12.8g

# Optional Color

White




Pearl White



Pink Purple



# Comparison of Old and New Versions

Appearance			
	Model	BK4870	BK-07B
EMS Electrode	6-pole electrode		Double electrode
Suction Power	✓		✓
Red Light	✓		✗
Battery Capacity	1400 mAh		700 mAh
Charging Method	Type-C		Dedicated adapter + USB cable
Anti-slip Design	✓		✗
Applicable body part	Body + Face		Body



# THANKS FOR YOUR ATTENTION

Contact for any Question

Mail: [cscare2021@gmail.com](mailto:cscare2021@gmail.com)

Phone : 516 653 0666

