



PL038 RF 5 in 1  
EMS Micro Current LED Light Facial  
Body Anti Aging Massage Device  
Operating Manual

Prohibition

Do not damage, destroy, forcibly bend, pull, twist, and clamp the plug cord, or place heavy objects on it.  
(Otherwise, it will cause malfunction, electric shock or fire.)

Please avoid overloading the socket.  
(Otherwise, it may cause heat or fire.)

Prohibition from disassembling

Do not disassemble, repair, or modify the device and base.  
(For repairs, please consult the retailer where you purchased it.)

Fire warning

Do not place this device in high-temperature places such as near a fire or under the scorching sun. In addition, please keep it away from heating appliances.  
(Otherwise, it will cause heat, rupture or fire.)

Do not put this device into fire and do not charge, use or place it in a heated or high temperature environment.  
(Otherwise, it will cause heat, rupture or fire.)

Be sure to confirm the shape and voltage of the plug before use.  
(Otherwise, it will cause heat, rupture or fire.)

Do not put this device in a microwave oven or high pressure container.  
(Otherwise, it will cause heat, rupture or fire.)

Compliance

Before use, be sure to take off precious metals, glasses and contact lenses.

When you use this device for the first time or have not adapted, be sure to start from the first speed.

If the device has not been used for a long time, please check whether it can start normally before use it.

When people who are undergoing dental treatment or have their teeth filled feel too irritating, please lower the speed and adjust the position of use.

Do not press nursing head hard on the skin, and do not stay in the same part continuously for more than 3s.

When caring for the skin around the eyes, close your eyes before use.

f your skin or body feels pain or abnormalities during use, please stop using it immediately and consult your doctor.

Part Identification

4

- In order to use this device safely and correctly, please read the Operation Manual carefully.
- After reading, please be sure to keep the Operation Manual properly for easy access at any time.
- When trans ferring this device to another person, be sure to attach the Operation Manual.

**Safety instructions**

- Please be sure to carefully read "safety Warnings" before use.
- Please use this device following the prescribed steps.
- Please stop using this device immediately if you feel uncomfortable.
- If use in a wrong me thod, it may cause adverse effects on your body.
- The other precautio ns specify im portant safety-related contents, please follow instruction.

Note

Danger

Improper use of this product will lead to serious casualties.

Warning

Improper use of this product may lead to serious casualties.

Attention

Improper use of this product may lead to personal injury or property damage.  
#Property damage includes the damage to the house, furniture, poultry and pets.

Sign

Refers to the prohibited matters.

Refers to the matters must be observed.

Refers to the matters must be noted.

Safety symbols

Type B applied part

Attention, see operating manual for use

Danger

use this device. Otherwise, accidents and physical discomfort may be caused due to the incorrect operation of medical electronic equipment.

Implantable medical electronic devices such as pacemakers.

Artificial heart-lung machine and other electronic equipment used to sustain life.

Wearable medical electronic equipment such as electrocardiograph and infusion set.

1

Function

Mode	Applicable period
BEAT -- intermediate frequency intermittent electric pulse	Suitable for weekly care
EMS -- Microcurrent muscle exercise	Suitable for daily care
TS--Low-frequency micro current	4-5 times a week
VR --Vibration Massage	Can be used alone or combined with other modes
LED -- 3 different lights 5 modes Red, green, orange, all lights on, all lights flashing.	Can be used alone or combined with other modes

■Timer function

During use, the device will emit a beep every 10 seconds.

■Automatic shutdown

The beauty device has the function of automatically powering off after 10 minutes.  
♥ After the machine works for 10 minutes, the device will emit a short beep 3 times, and turn off the power after a long beep.

Instructions

⚠ Please tear off the transparent protective film before using the product.

① Nursing guidance

1. Preparation before use

- This product should be used together with cosmetic liquid or cosmetic gel to avoid contact with products containing alcohol or acid substances.
- Please remove the metal objects (rings, necklaces, earrings, bracelets, etc.) from your body.
- Please remove makeup and cleanse before use to ensure that your skin is in good condition.
- When taking care around the eyes, remove the contact lenses first.
- When using in body care, it is recommended to take a shower before use to ensure that the skin is clean.

Removal of precious metal

Remove makeup and clean face

Bath to clean skin

5

⚠ Warning

- Do not use this device if you have a cardiac pacemaker, implanted defibrillator, or other implanted metallic or electronic device. Such use could cause electric shock, burns, electrical interference, or death.
- If you are in the care of a physician, consult with your physician before using this device;
- Do not use this device on your heart, head, mouth, pudendum or blemished skin areas;
- Do not use this device across your chest because the introduction of electrical current into the chest may cause rhythm disturbances to your heart, which could be lethal;
- Do not use this device over open wounds or rashes, or over swollen, red, infected, or inflamed areas or skin eruptions (e.g., phlebitis, thrombophlebitis, varicose veins);
- Do not use this device over, or in proximity to, cancerous lesions; - Do not use this device in the presence of electronic monitoring equipment (e.g., cardiac monitors, ECG alarms), which may not operate properly when the electrical stimulation device is in use;
- Do not use this device when in the bath or shower;- Do not use this device while sleeping;
- Do not use this device while driving, operating machinery, or during any activity in which electrical stimulation can put you at risk of injury;
- Do not use the device on children; Use this device only on normal, intact, clean, healthy skin;
- Do not use this device in high humidity areas such as a bathroom;
- Do not attempt to move the electrode heads while the device is operating; and
- Keep the device away from wet, high temperature and direct-sunlight place. Shady and dry place is good for device stock.

Prohibition

- Do not Overuse this device.  
(Please take care of this device in strict accordance with the use time and frequency of use specified in the Operation Manual. Otherwise, it may cause accidents or injuries.)
- Do not use this device for a long time on areas with thin skin or above bones.
- Do not use this device near flammable and explosive materials (Alcohol, thinner, gasoline, nail polish remover, spray, etc.)  
(Otherwise, it may cause explosion and fire.)
- Do not use this device with other cosmetic products.
- Do not use this device continuously at the same area, and be sure to while use it and moving it on the skin.  
(Otherwise, it may harm the skin or cause skin problems such as infection, inflammation, bleeding, or pigmentation.)
- Do not use this device for non-cosmetic purposes and in a damaged or malfunctioning state.

2

2. Charging

When using it for the first time (for the first time or not used for a long time), please make sure to charge it first. And confirm that it has been fully charged (Four indicators all on after fully charged.) before use.

• Charging time & times of use

Charging time	Times of use
About 2.5 hours	About 90 minutes

- ※ Please keep the charging environment at 25°C-35°C.
- ※ This is the operating time to be used immediately after full charging.
- ※ The service time will vary with the use mode, intensity and environment.
- ※ The charging time will be different due to the residual power of the charging battery and the use environment.

• Charging method

- 1). Insert the one socket of the charging line into the charging port of the device and the other socket into the USB port of the charger.
- 2). Plug the charger into the household power socket.

- 3). After charging is finished, pull out the charger from the socket, and then pull out the socket of charging line from the main unit.

- ※ This product only uses the provided charging line or the charging line and charger that meet the product specification requirements.
- ※ When using, plug the charger directly into the socket. No extension line is allowed to be used.

Please be sure to confirm the plug shape and voltage (charger output is DC5V= 1A) before use, otherwise it may cause body heating, fire or failure.

6

Forbidden To Use In A Humid Place

- Do not use or store this device in a bathroom or high humidity place
- Please avoid using this device at drinking places, bathroom, or while taking a bath.  
(If the device and charger are exposed by water, it may cause electric shock, short circuit or fire.)
- Try to avoid water on the device AC adapter and USB cable  
(Do not soak or wash it with water). (If the device and charger are exposed by water, it may cause electric shock, short circuit or fire.)
- Do not touch the device with wet hands.  
(Otherwise, it may cause electric shock or injuries.)

Uselimit

- Do not use this device near flammable materials, such as thinner, volatile oil, delustering water, spray.
- Do not use non-dedicated AC adapters.  
(Otherwise, it may cause malfunction, electric shock or fire.)
- Do not connect the AC adapter to other products.  
(Otherwise, it may cause malfunction, electric shock or fire.)
- Do not allow children to use this device, and avoid the disabled to use this device alone. Minors are requested to use this device with the consent and accompaniment of their guardian.  
(Otherwise, it may cause accidents or injuries.)
- Please do not allow people who cannot express their wishes and are very insensitive to skin perception to use this device  
(Otherwise, it may cause accidents and injuries.)
- Keep the device/AC adapter out of the reach of children. Do not allow children to use or play with this device  
(Otherwise, it may cause accidents or injuries.)
- Be sure to use 100-240V alternating current.  
(Otherwise, it may cause electric shock, short circuit or fire.)

Prohibition from dropping

- Do not place this device in an unstable place, do not drop it or bump the device and charger.  
(Otherwise, it may cause malfunctions or injuries.)
- Do not place heavy objects on the device or charger, and do not drop or bump this device during use.  
(Otherwise, it may cause injury, product failure, etc.)

Unplugging

When plugging or unplugging the AC adapter, be sure to hold the front plug for operation, and do not hold the cable for operation.  
• Do not use this device if the AC adapter is damaged or loosened after being plugged into the socket.  
• If there is an earthquake, thunder happen, or smoke from the equipment (host, base, AC adapter), abnormal heat (it is normal for this device to have a certain heat during use and charging), stink or deformation, please stop using this device, turn off the power, and unplug the AC adapter from the socket.

3

② Specific use

1. Start the beauty instrument

Press any buttons on the panel to turn on the device according to needs.

1) Start BEAT/EMS /TS mode

A. Power on, strength adjustment, power off	Turning on: after turning on, press "BEAT" or "EMS" or "TS" button to choose preferred mode with corresponding indicator on and beep warning tone, then you can start use. Level adjustment: press the function button (BEAT/EMS /TS) to adjust strength on certain mode, with 4 levels respectively.
B. Use with LED mode	After choosing specific function, adjust suitable strength level, and then LED lights. a. Press " " once, blue, orange and red light on; b. Press " " again, blue light on; c. Press " " again, orange light on; d. Press " " again, red light on; e. Press " " again, blue, orange, red light on; f. Press " " again, LED light off.  Red: for improving skin fatigue status caused by aging; Blue: for acne skin that is oily and with coarse pores; Orange: for improving skin tone and whitening; All light on: for comprehensive skin care; All light flashing: waking up skin before caring finished.
C. Use with VR mode	Under specific function, you can choose suitable strength level and LED light, and together with VR mode if you need. a. The first time you press the VR button, the VR mode is turned on. b. If you press the VR button a second time, the VR mode will be closed.

♥ BEAT/EMS/TS mode can be used with LED mode or VR mode or at the same time.

2) Turn on VR mode

Press the VR button. With a beep, the VR mode will start. If you press the VR button again, the VR function will be turned off.

♥ The VR mode can be used alone or in combination with any mode.

7

2.Skin care

1) Cycle care BEAT/TS

BEAT mode 5 minutes+ TS mode 5 minutes.

♥After cleansing, apply the gel. If the skin is dry in the middle of the operation, continue to apply the gel or beauty liquid.

♥Skin rejuvenation method: massage the skin in a circular motion along the skin texture.

1-2 times a week

①

②

③

Hold the beauty instrument and slide it from under the chin to behind the ear, then massage to the shoulder; then gently lift the neck upward.

Slide along the facial curve from the chin to the temples, from the nose to the temples, then from the corner of the eye to the end of the eye to the temples.

Slide your eyebrows to the hairline for a lift massage.

2) Daily care EMS

Once a day

Revitalize skin for 5-10 minutes every day.

♥After cleansing, apply gel. If the skin is dry in the middle of the operation, continue to apply gel or beauty liquid.

♥Skin rejuvenation method: massage the skin in a circular motion along the skin texture.

♥Place the beauty instrument in ● in sequence for no more than 3 seconds, and then slowly slide to the next part. Steps A to F are performed 2 to 3 times for each ● part.

②

②

③

②

②

3) Eye care TS

4-5 times a week

TS mode 5-10minutes.

♥After cleansing, apply the gel. If the skin is dry in the middle of the operation, continue to apply the gel or beauty liquid.

♥Nursing technique: Along the skin texture, the small massage head moves slowly against the skin.

①

Press the skin of the face with a small massage head to fully touch the eye lines, and slowly lift 1 cm vertically and then stay for about 3-5 seconds until the eye lines are stretched and tightened. Repeat 5 times and do the next one.

②

The small massage head gently massages the dark circles and bags under the eyes with moderate strength and slow movements, following the principle of moving from the inner corner of the eye to the outer corner of the eye. Repeat 8-10 times as appropriate.

③

Press the skin of the face with one massage head, fully contact the outer corner of the eye, and slowly lift up vertically, repeat 8-10 times as appropriate.

	<ul style="list-style-type: none"><li>• Do not use on Adam's apple, eye balls, mucosa and more positions.</li><li>• During using, do not stay in a certain position or press the skin hardly.</li></ul>
	<ul style="list-style-type: none"><li>• Repeat the operation on the left and right sides of the body and face several times.</li><li>• For the first time use, it is advised to start from the L level, and then choose proper strength levels according to your needs.</li></ul>

4) Body care

Body care EMS+LED+VR

Upper arm

The use time of one arm is about 10 minutes

1. Massage in circular motions under the breasts and under the armpits.

2. Lift and massage from bottom to top, repeating the massage on the left and right side several times.

Waist and abdomen

The total usage time is about 10 minutes

1. Massage in a circular motion around the navel.

2. Massage in circular motion from the waist down.

Hip

The total usage time is about 10 minutes

1. N-shaped massage of hips from inside to outside in a circular motion.

2. Lift from the bottom to the top, and massage left and right alternately several times.

Leg

The use time for one leg is about 10 minutes

1. Massage in circular motion from bottom to top.

2. First thighs, then calves, alternately massage left and right side several times.

Maintenance

- It is normal that the color of the massage head may change after the cosmetic liquid is improperly matched or cleaned improperly, or the beauty instrument is used for a long time.
- Before cleaning, please turn off the power and unplug it.
  - Protect the appliance from moisture.
  - This product is not waterproof. Do not immerse the appliance in water or rinse it with a faucet.
  - After use, wipe more than 3 times with dry paper towel with a small amount of water, or wipe with non-alcohol wet paper towel, wipe the beauty device and keep it in the dry place.
  - Do not use sharp or corrosive, volatile thinner cleaning products, otherwise it will cause discoloration or damage to the beauty instrument.

Storage

- When storing the product, avoid a hot or humid location (such as inside a car or bathroom) and keep it away from direct sunlight.
- When not in use, keep in a safe place to avoid tripping on it and getting injured.
- Keep out of the reach of small children and store away from pets.
- Do not store near fire.
- Before using the product after a long period of disuse, check and make sure there are no problems with the product.
- If it is not used for a long time (more than 2 months), please take the product out every two months and store it after it is fully charged.

Fault analysis and elimination

In case of any malfunction during the operation, please carry out following inspections. If failures cannot be eliminated or the product cannot operate properly, please contact Kingdom's after-sales service center. Do not disassemble and attempt to repair it by yourself.

Situation	Possible Reason	Measure
Charging light is not lit	Battery is dormant	The battery needs 6-30 minutes back to working. It is a normal phenomenon and needs to be fully charged before use.
	The charging wire or charger has been damaged.	The battery needs 6-30 minutes back to working. It is a normal phenomenon and needs to be fully charged before use.
Can't turn on the phone or turn it off after a "di" beep	Low battery	Connect the beauty instrument to the charger and fully charge it before using it.
During the use, there may be a feeling of needling/electric current/heating	It is not a fault but a normal effect of EMS or micro-current on skin. If you feel it too violently and cannot bear it, please use it at lower level or after fully moisturizing your skin with lotion before using it.	

Waste Treatment

The marking on this product or packaging material or accessory indicates that the material is made from a recyclable material. Separate waste into different recycling levels to facilitate the recycling of useful materials.

This appliance is equipped with lithium battery inside, so:

- Before the appliance is discarded, the battery must be taken out of it;
- When the battery is taken out, the appliance must be shut down and turned off;
- The battery shall be disposed of safely;
- The user cannot change the battery by himself. It is recommended that you send the product to the official collection point so that the professional removes the battery.

Remove the battery:

- When you use the tool to open the product and when you discard the battery, please take necessary safety measures.
- Take off the decorative cap with a tool;
  - Take out the four screws with a Philip's head screw;
  - Take off the upper cap after prizing it with a tool;
  - Unscrew the 2 screws from the PCB with a Philip's head screw and take the PCB out;
  - Unplug the battery lead and take off the battery.

Decorative cap

Specification

Product Name: RF 5 in 1 EMS Micro Current LED Light Facial Body Anti Aging Massage Device

Model:PL038

Rated Voltage:3.7V

Rated Power: 8W

Charging Power Supply: 5V 1A

Net Weight: about 134g

Product Dimensions: 50x49x222mm

Packing List: 1 main unit, 1 charger, 1 charging cable, 1 operating manual.

♥ Due to the improvement of the products, some specifications may be changed without prior notice.

Warranty exclusions:

This warranty does not apply to any costs, repairs, or services for the following:

- 1.Repairs necessitated by use other than normal home use.
- 2.Damaged resulting from misuse, abuse, accidents, alternation or improper installation.
- 3.Corrective work necessitated by repairs made by any-one other than Prospera authorized service technicians.

Non-warranty service is provided on a "per incident" basis. Consumer will verify that the product has failed and provide information for servicing the unit. All app-licable repairs, parts, shipping, handling, local tax and "per incident" fees will be charged for non-warranty repairs and support calls.

For warranty service, please call 925-292-5978, M-F:8:00am to 4:00pm, Pacific Standard Time. Or send email to sales@prosperacorp.com.

Consumer may also mail service request to this address: 392 West Larch Road # 31Tracy, CA 95304-1644

About Prospera

Prospera develops high quality, innovative, and affordable massage products, ranging anywhere from electronic pulse massagers to jade massagers. Prospera's products are available through leading department stores, consumer electronics stores, and major online retailers. Learn more about Prospera at www.prosperacorp.com

Prospera Corporation

392 West Larch Road # 31Tracy, CA 95304-1644

Tel:(925) 292 5978

Email:sales@prosperacorp.com



Warranty card



Warranty Registration Form

Product Name: \_\_\_\_\_

Product Model: \_\_\_\_\_ Purchase Price: \_\_\_\_\_

Date of Purchase: \_\_\_\_\_ Retailer: \_\_\_\_\_

First Name: \_\_\_\_\_ Last Name: \_\_\_\_\_

Address: \_\_\_\_\_

Cit: \_\_\_\_\_ State: \_\_\_\_\_ Zip: \_\_\_\_\_

Phone: \_\_\_\_\_ Email: \_\_\_\_\_

Signature and date: \_\_\_\_\_

Please mail to:  
392 West Larch Road # 31Tracy, CA 95304-1644  
Tel:(925) 292 5978  
Email:sales@prosperacorp.com