

Warranty Information

original purchaser of the product.
Promotion of this product is a co-

limited to the original consumer
date.

materials or workmanship.
shipment, unreasonable use,
authorized personnel or other
this warranty doesn't cover the
limited to; LCD Screens, glass
it is purchased and operated in
been used in violation of written
been altered or modified or, to
removed, altered, defaced or

expressed or implied and no
any other liability in connection
failure under any theory of tort,
gross negligence, strict liability,
Manufacturer's / Distributor's

period, a product with a defect
model (at manufacturer's option).
of the 90 day warranty period
such repair or replacement.
should first contact Nuvomed Inc
886, Monday to Friday 10am to

during incoming shipment. As
it, and insure it with the carrier.
me, return address and daytime
of your sales receipt or other
not be accepted.

NuvōMed®

Hydra-Gel Therapy Sleeve

CTSS-12/0209 | CTSL-12/0210 | CTSX-12/0211



Instruction Manual

Please read all instructions carefully and retain for future use

NuvōMed®

Getting Started

- Remove the product from the bag.
- Please read all instructions and retain all packaging until you have checked and counted all of the parts.

⚠ Caution

FOR EXTERNAL USE ONLY!

- Use only as directed.
- Always check temperature of product before use.
- Improper use may cause injury or burns.
- Discontinue to use and contact your physician if irritation occurs.
- Do not use this product on person with sensitive skin, decreased skin sensitivity, nerve damage or broken skin.
- Not recommended to use on face.
- Do not use while sleeping.
- Do not use if the product is damaged.
- Do not use if wet or if lining is damaged or gel exposed.
- Do not sit on, against, crush or destroy this product.
- Do not place it under your body.
- For Hot/Cold treatment, burns can occur, check skin under product frequently.
- Do not use on infants.
- Not to be used by unattended children or incapacitated person
- If any medical condition exists such as diabetes or if you are pregnant, consult a physician before use.

How to Clean

- To clean, use a damp cloth to mildly rub the outer layer of fabric.
- Recommended air dry after cleaning before storage.
- Lay the product out on a clean towel or flat surface to air dry.
- Choose an area away from direct sunlight to dry your product.

Nu

⚠
BAG
THE

How

Hot

⚠
MIC
DO NOT

Cold

⚠
DO NOT