

<p style="font-size: 24px; margin: 0;">105.00 mm</p>	<p style="font-size: 24px; margin: 0;">285.00 mm</p>	<div style="text-align: center;"> <p style="font-size: 24px; margin: 0;">LUNA'MOUR</p> <p style="font-size: 18px; margin: 0;">MicroMini Face Beauty Tool</p>  <p style="font-size: 18px; margin: 0;">User Manual</p> </div>	<p>Introduction Welcome to a new era of beauty with the LUNA'MOUR MicroMini Face Beauty Tool — a pocket-sized powerhouse designed to bring professional-level results into the palm of your hand.</p> <p>Crafted with elegance and precision, the MicroMini harnesses the power of advanced micro-current technology to gently awaken your facial muscles, restore youthful vitality, and encourage natural radiance. Each treatment feels like a mini workout for your skin — sculpting, lifting, and energizing from within — while remaining soothing, silent, and effortless.</p> <p>With every glide, the MicroMini works in harmony with your body's own bioelectricity, helping to firm contours, smooth fine lines, and enhance your natural glow. Its compact design makes it the ultimate companion — whether at home or on-the-go.</p> <p>Key Features</p> <ul style="list-style-type: none"> • Micro-Current Technology: Low-level currents mimic the body's natural bioelectricity to stimulate collagen and elastin. • Facial Exercise Training: Tones and strengthens facial muscles for improved definition. • Compact & Lightweight: Just 35 grams, perfect for travel and daily use. • Silent Operation: Less than 50 decibels for a relaxing session. • Rechargeable Battery: Built-in 3.7V battery (170mA). • Precision Design: Ergonomic dual-contact points for effective conductivity. <p>How to Use</p> <ul style="list-style-type: none"> • Press the power button once to start at Level 1. Press again to switch levels. • To turn Vibration Mode on, double-click the button to activate the vibration feature. • To turn off the device, press and hold for 2 seconds. • Charging: Red light ON while charging, turns OFF when fully charged. <p>Tips for Best Results</p> <ul style="list-style-type: none"> • Cleanse and dry face before use. • Apply water-based conductive gel or serum. • Glide device gently along jawline, cheeks, and forehead. • Use for 5–10 minutes per session, 3–5 times per week. • Consistency is key for best results. 	<p>Isafety & Precautions</p> <ul style="list-style-type: none"> • Do not use on broken, irritated, or infected skin. • Avoid if you have a pacemaker, metal implants, or epilepsy. • Not recommended during pregnancy without medical approval. • Keep out of reach of children. • Stop use if you experience discomfort. <p>Maintenance & Care</p> <ul style="list-style-type: none"> • Wipe device with a soft cloth after each use. • Do not submerge in water. • Store in a cool, dry place. • Charge fully before first use. <p>Product Specifications Model Name: LUNA'MOUR MicroMini Face Beauty Tool Material: ABS Plastic Dimensions: 57 x 52 x 24 mm Weight: 35 g Battery: 3.7V / 170mA rechargeable lithium battery Frequency: 185 Hz Noise Level: < 50 dB Patent Number: 202030103955.6 UPC: 850064305104 SKU: LL-MM</p> <p>Warranty & Client Care This device is covered by a limited 1-year warranty. For warranty services, please contact our Client Care team at hello@lunamour.co. Warranty covers only manufacturing defects when purchased through authorized partners.</p> <p>Disclaimer This product is not intended to diagnose, cure or prevent any disease. LUNA'MOUR makes no claims regarding the ability of this product to treat medical conditions. Always consult with a qualified health professional for any concerns.</p>
--	--	--	---	---

<p style="text-align: center; font-weight: bold;">Warnings:</p> <p>CONTRAINDICATIONS: Ask your doctor prior to using if you: Are taking prescription or over the counter drugs known to cause light sensitivity. Are, or suspect you are, pregnant. Have any suspicious or cancerous lesions.</p> <p>Please read the instructions carefully before use and use this product in the correct order.</p> <p>Please keep out of reach of children. This device is designed for use on adults over 18 years of age and not recommended for use on children under the age of 12. Please do not use if you are experiencing the following: Heart disease, high fever, infectious diseases, malignant tumors, skin wounds, acute diseases, using medical appliances, skin diseases, pregnant women and/or breastfeeding, excessive sunburn, particularly sensitive skin, and those who are in the hospital.</p> <p>Please consult your doctor prior to use if you are, or suspect you are, pregnant. Please consult your doctor prior to use if you have, or suspect you have, cancer. Do not use device on open wounds. Do not use this device on skin that is painful to touch, raised (higher than the surrounding skin), crusting, inflamed, or bleeding. This may be a sign of a more serious condition and should be treated by your doctor. No modification of this device is allowed.</p> <ul style="list-style-type: none"> -Do not place magnets, metal or conductive objects on the body of the device as it may cause malfunction. -Do not rinse the instrument with water, use protective cloth provided to wipe clean after use. -Do not wipe with alkaline detergent, alcohol, thinner, gasoline, nail polish remover, etc. -If not in use for a long time, please make sure that the instrument is fully charged. -Please stop using it immediately when there is abnormal heat, so as not to cause burns and fires and other safety issues. -When the power cord is damaged, please do not use it to avoid short circuit and fire -Do not place near flammable materials (spray, oil, natural gas, etc.) -DO NOT submerge the device into water as it will malfunction and/or cause injury. -If product has abnormal noise, odor, or smoke, do not disassemble it at will, please contact the customer service email: hello@lunamour.co 	<p style="text-align: center; font-weight: bold;">Important Safety Information:</p> <p>-If the product is faulty, please contact the customer service center, please do not disassemble or modify the machine by yourself. Do not operate device where it would come in contact with water. Do not use or charge device with power adapter that is damaged and/or faulty. Do not harshly bend or abuse the wire of USB charging cord. Do not attempt to take the device apart. Device is not serviceable by user. This will void any warranty.</p> <p style="text-align: center;">Stop using the device immediately if you think it is defective. Please contact Client Care at hello@lunamour.co</p> <p>Product Representation:</p> <p>LUNA' MOUR is constantly improving its products with the latest technology. As a result, the actual product may differ slightly from the product shown in the images on the packaging and in the user guide.</p> <p>Skin Patch Test: Prior to using the device for the first time, it is recommended that you conduct the following skin test on your forearm:</p> <ol style="list-style-type: none"> 1) Keep your forearm dry and clean. 2) Turn on the device which will activate level one of the microcurrent mode, and conduct a three-minute treatment. Gently press device over your skin so it is in full contact. If rash, itch or irritation occurs on the tested skin, please stop the test immediately. If skin feels abnormal, also stop use immediately. 3) If the area is not red and no negative condition or fever appears after 24 hours, then you may keep using the device. 	<p style="text-align: center; font-weight: bold;">MicroMini</p> <p style="text-align: center;">Client Care: For information regarding device usage, warranty, and additional information, please contact Client Care at hello@lunamour.co</p> <p style="text-align: center;">Disclaimer: This product is not intended to diagnose, cure or prevent any disease. LUNA' MOUR makes no claims, representations, or warranties regarding the ability of this product to cure any mental, physical, or skin conditions. A qualified health professional should be consulted with regards to any condition requiring medical attention.</p> <p style="text-align: center;">www.lunamour.co hello@lunamour.co</p> <p style="text-align: center;">Item Number : S022 LL-MM</p>
---	---	---