



fat Freezer

Sculpt your body by freezing your fat cells

INSTRUCTION MANUAL



Dear Valued Customer,
Congratulations on your Fat Freezer purchase. We are certain that it will become an essential part of your toning routine. The Fat Freezer speeds up the natural elimination of fat cells and is a machine, intended for use by adults, desiring to decrease subcutaneous fat deposits in specific areas of the body.

Warning:
Don't be alarmed if during use, the Fat Freezer machine may heat up as this is a normal process to allow the machine to ventilate after extended use.

Thank you

Beautyko LLC. • Phone: (516) 653-0666 • Email: onlineorders2017@gmail.com • www.beautyko.com

Features:

The FatFreezer uses scientifically proven cold lipolysis technology that works by delivering energy (cold) through an area of the skin that helps dissolve the fat cell deposits beneath that area by "freezing" them away.

FatFreezer provides a non-invasive safe and simple way of reducing pockets of fatty tissue in designated areas of the body.

Important safety Instructions:

WARNING! No modification of this equipment is allowed.

FatFreezer is powered through a power adapter which converts AC power from a wall outlet to the DC power required by the FatFreezer.

Your FatFreezer comes with a universal AC-DC adapter. That means that you may connect the power cord to any 100V-200V as well as 220V-240V outlets without setting switches or using power converters.

Defective FatFreezer :

Stop using the device immediately if the device does not work properly. Please attach a note with the following "Do not use this FatFreezer" to prevent any further use. This defective device should be stored in a safe place until it has been returned for replacement.

Environmental conditions:

Transport & storage

Ambient temperature: -20°C to 60°C (-4°F to 140°F)

Relative humidity: 0% to 85%

Atmospheric pressure: 700hPa to 1060hPa

Normal use

Ambient temperature: 20°C to 25°C (68 to 77°F)

Relative humidity: 30% to 85%

Atmospheric pressure: 700hPa to 1060hPa

Electrical caution:

DO NOT connect the AC power cord to an AC outlet before connecting the DC plug to the fatFreezer because doing this may damage the AC-DC adapter.

IMPORTANT! Damage may occur if you use a different adapter to power the FatFreezer or if you use FatFreezer adapter to power other electrical devices. If there is smoke, burning scent, or extreme heat coming from the AC-DC adapter. You may damage your fatFreezer with a faulty AC-DC adapter.

WARNING! The power adapter may become warm to hot when in use. Be sure not to cover the adapter, and to keep it away from your body.

IMPORTANT! Unplug the power adapter or switch off the AC outlet to minimize the power consumption when the FatFreezer is not in use.

Health warnings:

1. WARNING! Do not use this product if you are pregnant, a child under the age of 18 years of age, incapacitated or immobile person or anyone who cannot clearly understand the instructions and/or operate the controls.
2. WARNING! Do not use this product if you have any heart-related illnesses, suffer from epilepsy, fever, inflammation, dermatitis (acne, eczema, herpes), warts open or recently stitched wounds, serious or progressive illnesses (cancer, AIDS, hepatitis, blood diseases, etc.), blood coagulation problems or mental illness.
3. WARNING! Do not use the FatFreezer on your head, eyes, neck (throat), breasts, genital area, and inside your mouth.
4. WARNING! Do not use the FatFreezer without protective pad.

Possible side effects:

Temporary but rare side effects include bruising, numbness, redness and irritation of the skin where the fatFreezer has been placed. These side effects are rare and only temporary.

If skin irritation occurs after treatment, please discontinue use immediately. IF symptoms persist contact a physician.

Proper use guidelines:

WARNING! It is unsafe to use the FatFreezer, its accessories, detachable parts and materials in ways not described in the instructions for use.

Before use:

Make sure that the FatFreezer works correctly by turning on the machine briefly before using it.

Make sure that the FatFreezer is properly connected.

Make sure that the all accessory parts are properly attached.

Check and maintain the cleanliness of the parts in direct contact with your skin via the protective membrane.

Avoid placing the device near a strong magnetic field or static electricity.

Avoid placing the device near high frequency equipment.

During Use:

WARNING! Do not clean or maintenance the FatFreezer while the equipment is in use.

WARNING! Do not position the FatFreezer so that it is difficult to disconnect the power cord of the adapter. Maintain the proper positioning of the FatFreezer with the CONTROL PANEL face upwards. This will allow quick access to the OFF button and quick access to removing the DC plug from the device.

Select the desirable "Time" option. Once the time has been reached, turn off the FatFreezer and disconnect it from the power cord.

NEVER use the FatFreezer twice (back-to-back) on the same area of your body in a single session.

WARNING! Stop using the FatFreezer if you experience unusual pain. Immediately press the "OFF" button to stop the machine and disconnect it from the power cord.

WARNING! Stop using the device if you notice any abnormalities (for example: liquid inside the device, inability to turn on as intended, or an unusual sound with the fan).

After use:

Let the FatFreezer set for 10 minutes after turning it off and disconnecting it. This allows time for the metal plate to return to normal (room) temperature.

Clean the FatFreezer by wiping over it with a soft cloth.

You may use an organic solvent (for example: antiseptic solution, etc.) to clean the device.

After using the FatFreezer, you MUST wait 2 hours before using it again.

Keep the original box and wrapping to store the FatFreezer and to transport it

Parts Identification:

- | | | |
|-------------------------------|-----------------------|-----------------------|
| 1. control panel | 2. 30-min time button | 3. 60-min time button |
| 4. "OFF" button | 5. Belt loop handle | 6. Belt loop handle |
| 7. Metal plate (applied part) | 8. Cooling fan | 9. DC connection |

Accessories

- | | | |
|--------------------------|------------------------------|--------------------|
| 1. Velcro belt | 2. AC-DC adapter/ power cord | 3. Protective pads |
| 4. Recommended diet plan | 5. Fat measuring device | |

Setup instructions:

1. Loop one end of the Velcro belt around the handle of one end of the FatFreezer machine.
2. Apply a protective pad on the desired area of your body.
3. Carefully place the FatFreezer machine (metal plate) over the protective membrane.
4. Wrap the Velcro belt firmly around the body part and apply the other end around the opposite loop attach the FatFreezer firmly.

WARNING! To avoid strangulation caused by cutting off or impeding your blood circulation. DO NOT tightened the FatFreezer excessively tight around your body part. The Velcro belt is used to hold the FatFreezer firmly to avoid slipping, and NOT to obstruct blood circulation.

WARNING! If you feel numbness on the body part where the FlimFreezer immediately detached the device from your body

part.

5. Plug the AC/DC adapter/ power cord into the fatFreezer and then into an electrical outlet.

Operating instructions:

6. Press either the '30' or '60' time allotment button to turn on the FatFreezer machine. (Note: 30 minutes for toning area/60 minutes to treat fat cells in that area.)

Button Function

OFF: Pressing this button turns off the FatFreezer. Press this button when you want to interrupt timed operation (30 or 60 minutes).

30: Pressing this button enables the FatFreezer to operate for 30 minutes. It will automatically shut off once the 30 minute is reached.

60: Pressing this button enables the FatFreezer to operate for 60 minutes. It will automatically shut off once the 60 minute is reached.

IMPORTANT! A green light turns on after pressing the '30' or '60' button to confirm that the FatFreezer machine is working properly.

7. When the designated time is reached, the FatFreezer will beep and shut off.

8. Unplug the FatFreezer from the AC outlet.

9. Let the FatFreezer sit for 10 minutes to return to normal temperature before cleaning and storing it away.

10. Gently massage the area of the body, using circular motions for at least 2 minutes.

WARNING! You must wait 4 weeks before applying the FatFreezer to the same area.

WARNING! If skin irritation occurs after treatment, discontinue use immediately. If symptoms persist, contact a physician.

WARNING! Cleaning procedures:

Dirt, grime and dust can accumulate and negatively affect the performance of FatFreezer and decrease service expectancy.

Use a soft cloth to wipe down the FatFreezer after each use. The entire device can be cleaned.

Periodically, you may use a dust-removal blower to clear the fan area of dust accumulation.

Repair :

If your FatFreezer malfunctions during the warranty. DO NOT attempt to open the device, as this may damage the delicate electrical components in the device.

Disposable parts

The protective membranes are the only disposable accessories used with the FatFreezer.

Protective membranes should be discarded after use. The protective membranes can be discarded in your household receptacles/ waste containers in the same manner you would discard normal waste.

Replaceable parts

There are no replaceable arts that require installation by the consumer to repair the FatFreezer. The manufacture will replace a machine that has been assessed as defective.

Storage precautions

WARNING! Do not leave the device where it could be influenced by vibration or shock.

WARNING! Do not leave the device in a dusty, salty or sulfuric environment. These elements can negative affect the performance of FatFreezer and decrease service expectancy.

Troubleshooting guidelines:

Does Not Turn On

- 1.FatFreezer is properly connected to electrical outlet. Power cord is properly connected to the AC-DC adapter, and that adapter is connected to the FatFreezer, and power cord is connected to electrical outlet.

- 2.Ensure that the electrical outlet has power.

Does Not Turn Off

- 1.If you press the "OFF" button to stop the FatFreezer and it does not turn off, immediately disconnect the AC-DC adapter from the machine.

- 2.You can then disconnect the power cord from the electrical outlet.

Beautyko LLC
Phone: (516) 441-5299
customerservice@beautyko.co
www.beautyko.com



BK1136