

UNDER EYE REJUVENATOR

Item No. 208311

User Guide



Thank you for purchasing the Under Eye Rejuvenator. Please take a moment to read this guide and store it for future reference.

INTRODUCTION

The Under Eye Rejuvenator combines vibration, heat and electrical muscle stimulation [E.M.S.] to soothe and revitalize the areas under your eyes. It promotes circulation to relieve fatigue and helps reduce the appearance of fine lines, wrinkles and dark circles.

FEATURES




- Comfortable, ergonomic frame
- Cordless and portable
- Vibration mode with 3 speeds
- Hot compress mode
- Infrared LED therapy
- E.M.S. pulse therapy with adjustable intensity
- Charges via USB

CHARGING THE DEVICE

- This device arrives partially charged. However, for best performance, please charge the device fully prior to first use.
- Plug the included cable into the USB charging port.
- Connect the USB charging cable to a computer, USB power bank or AC adapter [USB to AC adapter not included].
- The charging time is approximately 2 to 3 hours.
- While charging, the green LED will be flashing. When it is fully charged, the LED will turn solid green.

DESCRIPTION OF PARTS

- A. Infrared therapy lamps
- B. HEAT button
- C. Power button
- D. E.M.S. Pulse button
- E. Stainless steel massage surfaces
- F. USB charging port



GETTING STARTED

1. Please charge the device fully before first use.
2. Wash the face thoroughly before using this device.
3. After washing the face, apply a hot towel to the eye area before using this device.
4. Apply your favorite moisturizer to the under-eye area.
5. Remove the protective film covers from the stainless steel massage plates.
6. Place the included conductivity pads onto the stainless steel massage plates.
7. Place the device on your face so it rests comfortably and naturally.

VIBRATION MASSAGE MODE

Press and hold the **POWER** button for 2-3 seconds. The green LED will illuminate, and you will hear one beep, indicating that the device is in Vibration Massage Mode, level 1 [low speed].

Press the **POWER** button a **SECOND** time. You will hear 2 beeps, indicating that the device is in Vibration Massage Mode, level 2 [medium speed].

Press the **POWER** button a **THIRD** time. You will hear 3 beeps, indicating that the device is in Vibration Massage Mode, level 3 [high speed].

Press the **POWER** button a **FOURTH** time. You will hear 1 long beep. The vibration massage will stop, however the green LED will remain illuminated, indicating that the device is in standby mode.

To turn the device **OFF** completely, press and hold the **POWER** button for 2-3 seconds. The device will beep once and the green LED will turn off.

HEAT / INFRARED LIGHT THERAPY

While the device is in standby mode, press the **HEAT** button. You will hear 1 beep, indicating that the hot compress function is active.

Press the **HEAT** button a **SECOND** time. You will hear 2 beeps, and the red LED will illuminate, indicating that the hot compress **AND** the infrared light functions are **BOTH** active.

Press the **HEAT** button a **THIRD** time. You will hear 3 beeps, and the red LED will remain illuminated. In this setting, the hot compress will be turned off, but the infrared light therapy will remain on.

Press the **HEAT** button a **FOURTH** time. You will hear 1 long beep. The Heat / Infrared Light Therapy will stop, however the green LED will remain illuminated, indicating that the device is in standby mode.

E.M.S. PULSE MODE

Press and hold the **POWER** button for 2-3 seconds. The green LED will illuminate, and you will hear one beep, indicating that the device is in vibration massage mode 1 [low speed]. At this point, press and hold the **PULSE** button. You will hear one beep. This indicates that you have switched from vibration mode to E.M.S. Pulse mode. The default is level 1 [LOW intensity].

Press the **PULSE** button a Second time, you will hear two beeps. This indicates that the E.M.S. Pulse mode is on level 2 [MEDIUM intensity].

Press the **PULSE** button a third time, you will hear three beeps. This indicates that the E.M.S. Pulse mode is on level 3 [HIGH intensity].

Press the **PULSE** button a **FOURTH** time to cycle back to level 1. Continue pressing to cycle through the other levels.

To turn the device **OFF** completely, press and hold the **POWER** button for 2-3 seconds. The device will beep once and the green LED will turn off.

NOTE: *The E.M.S. Pulse Mode delivers a gentle electric microcurrent to the under-eye area, which produces a tingling sensation in the face. Always be sure to start your treatment at the **LOWEST** intensity until you become accustomed to E.M.S. Pulse Mode treatment. As you become more familiar with the sensation, you may wish to increase the intensity to Level 2 or Level 3. Moisturizing skin will improve the conductivity of the electric currents. Therefore, it is recommended to use the included conductivity pads or moisten your skin with lotion before activating the E.M.S. Pulse function.*

COMBO TREATMENTS

As you become more familiar with the Under Eye Rejuvenator you may wish to customize your routine by combining the various treatment modes:

Vibration + Heat / Infrared Therapy

Press the **POWER** button once to activate Vibration Massage Mode. While the vibration is activated, Press the **HEAT** button. At this point, the massage and heat functions are now operating simultaneously. You can adjust each function separately by pressing the **POWER** and/or **HEAT** buttons as described above.

E.M.S. Pulse + Heat / Infrared Therapy

Activate the E.M.S. Pulse Mode as described above. While the Pulse function is on, press the **HEAT** button. At this point, the E.M.S. Pulse and heat functions are operating simultaneously. You can adjust each function separately by pressing the E.M.S. Pulse and/or HEAT buttons as described above.

PLEASE NOTE: *It is not possible to combine Vibration Mode with E.M.S. Pulse Mode.*

RECOMMENDED USAGE

The recommended usage is 3 times per week.

SAFETY WARNINGS

- Not for use in excessively hot or dusty environments.
- Not for use in excessively wet or humid environments.
- This device is not a toy. Keep out of reach of children at all times.
- This device is for personal use only, and not intended for commercial use.
- Do not attempt to open or modify this device. There are no user-serviceable parts.
- Do not use this device while driving or operating machinery.
- Do not use this device for anything other than its intended purpose.
- Be sure this device is powered off and unplugged when not in use.
- Use only the charging cable that came with this device.
- When not in use, store this device in its original carton in a cool, dry location.
- Consult a physician before using any at-home massage device.

NOTE: *Some users may experience redness in the face after using this device. If this occurs, apply a warm towel to the face. The redness should subside after 30 minutes.*

THIS ITEM IS NOT RECOMMENDED for those who have facial trauma, severe skin inflammation, facial fractures, high/abnormal blood pressure or eardrum inflammation; those who are recovering from surgery or cosmetic surgery, sensitive to metal, under 12 years old or pregnant; and those who wear hearing aids.