

MILEX INSTA-FOLD FOLDABLE BED INSTRUCTION MANUAL:

- 1) The Foldable Bed when fully Stowed is only 8 inch in width
- 2) To Release the latch and press down on the middle joint position to expand the mattress
- 3) Once fully expanded secure the Bed and confirm all positions are locked in place.
- 4) To release, pull the latch underneath and easily fold together

① *Stowed state,
width of only 8inch*



② *Release the latch and
press down on the middle
joint position*



③ *The bed surface is fully expanded,
flat and comfortable.*

